CARESS FLOW

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
LO STUDIO CLINICO PILOTA

EUROPEAN REVIEW FOR MEDICAL AND PHARMACOLOGICAL SCIENCES

2018; 22: 8480-8486

Vaginal natural oxygenation device (VNOD) for concomitant administration of hyaluronic acid and topical hyperbaric oxygen to treat vulvo-vaginal atrophy: a pilot study

L. CONDEMI¹, J. DI GIUSEPPE², G. DELLI CARPINI², F. GAROIA³, A. FREGA⁴, A. CIAVATTINI²

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⁴Department of Surgical and Medical Sciences and Translational Medicine, Sapienza University of Rome, S. Andrea Hospital, Rome, Italy

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Abstract. - OBJECTIVE: This is a pilot study to evaluate the effectiveness of concomitant administration of hyaluronic acid and topical hyperbaric oxygen therapy (THOT) by a specifically designed medical device (vaginal natural oxygenation device, VNOD) in improving the symptomatology of postmenopausal patients with vulvovaginal atrophy (VVA).

PATIENTS AND METHODS: Women with diagnosis of severe VVA from September 2017 to May 2018 were included. Five biweekly administration of THOT and concomitant of hyaluronic acid were performed with a specifically designed medical device. In each occasion, the intensity of patient's symptoms (wellbeing such as absence of dyspareunia, vaginal dryness, vulvar and/or vaginal itching; vaginal burning; presence of fluid) was determined with a graduated scale from 1 to 6 and the vaginal elasticity and the vaginal wall epithelium appearance were also determined with a graduated scale from 1 to 5. The change in all parameters from baseline to end of therapy was evaluated.

RESULTS: Twentyfive patients were considered for the final analysis. A significant improvement in well-being (0.3 vs. 5.1, $p < 0.001$), vaginal burning (0.2 vs. 5.1, $p < 0.001$), presence of fluid (0.6 vs. 4.9, $p < 0.001$), vaginal epithelium appearance (1.8 vs. 4.7, $p < 0.001$), and vaginal elasticity (1.1 vs. 3.8, $p < 0.001$) was observed between the first and the last therapy session. All the patients reported a recovery of their sexuality at the end of the five treatment sessions.

CONCLUSIONS: In this pilot study, the use of VNOD seems to be a valid treatment of VVA, resulting in a completely natural type of therapy well accepted by patients with immediate therapeutic effects and without side effects; these findings must be confirmed in a well-designed randomized controlled trial.

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CARESS FLUX
The most advanced in motion chair

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